

## TS 208 – TDA REAGENT (10 ml/vl)

### INTENDED USE

For determination of tryptophan-deaminase activity.

### COMPOSITION

Ingredients	Concentration
Ferric chloride	10.000gm
Distilled water	100.000ml

### Applied in Medium













TM 1639 - CHROMOGENIC UTI AGAR, MODIFIED

### INSTRUCTION FOR USE

Add 1-2 drops of TDA Reagent directly on suspected colony from TM 1639 - CHROMOGENIC UTI AGAR, MODIFIED. Observe for appearance of dark brown colouration around the colony within 10 – 30 seconds for confirming positive reaction.

### STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

**\*For Lab Use Only**  
Revision: 9 March., 2022