

TS 059 – DIPHTHERIA VIRULENCE SUPPLEMENT (Part A&B) (Double Pack) KIT

INTENDED USE

For isolation and presumptive identification of *Corynebacterium diphtheriae*.

COMPOSITION

Ingredients	Concentration
Part A	
Horse serum	100ml
Part B	
Potassium tellurite	1ml

(per Kit sufficient for 1000 ml medium)

Applied in Medium













TM 450 - TINSDALE AGAR BASE

INSTRUCTION FOR USE

Warm up the refrigerated contents of Part B and aseptically add 29 ml sterile distilled water. Mix thoroughly. Aseptically add warmed up (to 50°C) contents of Part A (Horse serum) and Part B (Potassium tellurite) to sterile, molten, cooled (45-50°C) TM 450 - TINSDALE AGAR BASE. For 10 ml of TM 450: 1.0 ml of Part A and 0.3 ml of Part B, is recommended. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative	 CE European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 10 Aug., 2024