

TMK 362S– BHI SUPPLEMENTED W/ 0.05% SPS

INTENDED USE

For detection of microorganisms associated with blood culture.

PRODUCT SUMMARY AND EXPLANATION

BHI Supplemented W/ 0.05% SPS is used as a blood culturing system for cultivating a wide variety of microorganisms since it is a highly nutritive medium. This culture medium complies with the specifications given by EN ISO 6888 and APHA. BHI Broth is a modification of the original formulation of Rosenow, where he added pieces of brain tissues to dextrose broth. It is an enriched non-selective broth medium that is useful in the cultivation of fastidious and non-fastidious microorganisms. This medium will also support the growth of aerobic microorganisms from a variety of clinical and non-clinical specimens. For selective isolation of fungi, addition of gentamicin and/or chloramphenicol is recommended.

COMPOSITION

| Ingredients | Gms / Ltr |
|----------------------------------|-----------|
| Calf Brain infusions from 200 gm | 12.500 |
| Peptone | 10.000 |
| Beef heart Infusion from 250 gms | 5.000 |
| Sodium chloride | 5.000 |
| Disodium hydrogen phosphate | 2.500 |
| Dextrose | 2.000 |
| SPS | 0.500 |

PRINCIPLE

The mixture of brain and heart infusions provides organic nitrogen, carbon, and vitamins. Dextrose is the carbohydrate source. A low concentration of dextrose is used to stimulate early growth. Sodium chloride maintains the osmotic environment. Disodium phosphate is the buffering agent in this medium and also helps neutralize the acids produced from the utilization of dextrose, thus maintain viability.

Sodium polyanethol sulfonate (SPS) is an anticoagulant and a surface-active agent which is widely employed as an additive to fluid blood culture media. It is generally considered to enhance the rate and speed of bacterial isolations by counter-acting the bacterial inhibitors of human blood. SPS is known to neutralize the bactericidal activity of fresh human serum and to inhibit phagocytosis.

INSTRUCTION FOR USE

1. Remove the plastic cap and disinfect the part of the rubber stopper which is now exposed.
2. Draw patient's blood with the sterile needle and syringe and transfer the blood sample immediately into the culture bottle by puncturing the rubber stopper with the needle and injecting the blood.
3. Venting may be required for aerobic culture and not in case of anaerobic cultures.
4. Incubate at 35-37°C for 18 -48 hours and further for 7 days to confirm negative results.

Note: BHI Supplemented W/ 0.05% SPS is a ready to use liquid media in glass bottle. The medium is pre-sterilized, hence sterilization is not required.



QUALITY CONTROL SPECIFICATIONS

| | | |
|---------------------------------|---|---|
| Appearance of the medium | : | Light amber colour, clear solution. |
| Quantity of Medium | : | 25ml / 50ml of the medium in glass bottle |
| pH (at 25°C) | : | 7.4 ± 0.2 |
| Sterility Check | : | Passes release criteria |

INTERPRETATION

Cultural characteristics observed after incubation.

| Microorganism | ATCC | Inoculum (CFU/ml) | Growth | Incubation Temperature | Incubation Period |
|--------------------------------|-------|-------------------|-----------|------------------------|-------------------|
| <i>Escherichia coli</i> | 25922 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |
| <i>Staphylococcus aureus</i> | 25923 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |
| <i>Staphylococcus aureus</i> | 6538 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |
| <i>Streptococcus pneumonia</i> | 6303 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |
| <i>Streptococcus pyogenes</i> | 19615 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |
| <i>Enterococcus faecalis</i> | 29212 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |
| <i>Neisseria meningitidis</i> | 13090 | 50-100 | Luxuriant | 35-37°C | 18-48 hours |

PACKAGING:

Aluminium capped bottles containing 25ml (Paediatric) or 50 ml (Adult) media.

STORAGE

On receipt, store bottles in the dark at 10 to 25° C. Avoid freezing and overheating. The medium may be used up to the expiration date and incubated for the recommended incubation times. Bottles from unopened packages can be used up to the expiration date. Opened bottles must be used immediately.

Product Deterioration: Do not use bottles if they show evidence of microbial contamination, discoloration, or any other signs of deterioration.

DISPOSAL

After use, prepared media, specimen/sample containers and other contaminated materials must be sterilized before discarding.

REFERENCES

1. Rosenow, Dental Research, 1, 205. (1919).
2. Rosenburg T. et al, J. Inf. Dis., 74, 131. (1944).
3. Mc Faddin J.F., Media for Isolation-Cultivation-Identification-Maintenance of medical Bacteria, Vol. I, Williams and Wilkins, Baltimore (1985).
4. Lennette, Balows, Housler and Shadomy (Eds.), Manual of Clinical Microbiology, 4th ed., ASM, Washington, D.C. (1985).
5. Ajello L., George L., Kaplan W. and Kaufman L., CDC Laboratory Manual of Medical Mycology, Atlanta, Ga.: US. DHEW, Center for Disease Control. (1966).
6. McDonough E., Geoge L., Ajello L. and Brinkman S., Mycopathol. Mycol. Appl.; 13, 113. (1960).
7. Selection for vancomycin resistance in clinical isolates of Staphy. haemolyticus: : R.S. Schwalbe, et al., J. Infect. Dis. 161, 45. (1990).



QTY.
Quantity

**LOT/
B. NO.**
Lot / Batch Number


Temperature Unit


Manufacturer


Best Before

GMP
Certification of
Good Manufacturing Practices

REF
Catalogue No.

EC REP MedMar GmbH
Buckenhorn 10,
49163 Mander, Germany
Authorized Representative


European Conformity


Consults Instructions for use :


QR
Code

IVD
For In Vitro Diagnostic Use

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 16th Feb. 2022

