

TS 020 – NALIDIXIC SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of *Pseudomonas aeruginosa* from clinical specimens.

COMPOSITION

Ingredients	Concentration
Nalidixic acid	15mg

(per vial sufficient for 1000 ml medium)

Applied in Medium

TM 2035 - CETRIMIDE AGAR BASE (W 1.3% AGAR)

TM 060 - CETRIMIDE AGAR BASE


TMV 060 - CETRIMIDE AGAR BASE (VEG.)

INSTRUCTION FOR USE

Rehydrate the content of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 1000 ml of sterile, molten, cooled (45-50°C) TM 2035 - CETRIMIDE AGAR BASE (W 1.3% AGAR) / TM 060 - CETRIMIDE AGAR BASE / TMV 060 - CETRIMIDE AGAR BASE (VEG.). Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 4 March., 2022