

TS 030 – UREA 40% (5 ml/vl)

INTENDED USE

For detection of urease activity.

COMPOSITION

Ingredients	Concentration
Urea	2g
Distilled water	5ml

(per vial sufficient for 100 ml medium)

Applied in Medium













TM 477 - DIAGNOSTIC STUART'S UREA BROTH BASE (UREA BROTH BASE)
TMV 477 - DIAGNOSTIC STUART'S UREA BROTH BASE (UREA BROTH BASE) (VEG.)
TM 1851 - UREA BROTH BASE (CHRISTENSEN)
TM 1196 - HEMMES MEDIUM
TM 1554 - KOHN TWO TUBE MEDIUM NO. 1
TM 925 - YERSINIA IDENTIFICATION BROTH BASE
TM 781 - MIU MEDIUM BASE
TMV 781 - MIU MEDIUM BASE (VEG.)
TM 394 - UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE)
TMV 394 - UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE) (VEG.)
TM 1113 - UREA AGAR BASE (CHRISTENSEN)
TM 2410 - UREA AGAR BASE (CHRISTENSEN) (IS:5887 (Part-I)-1976)

INSTRUCTION FOR USE

Warm up the refrigerated Urea Solution to room temperature and aseptically add 5 ml in 95 ml sterile, molten, cooled (45-50°C) TM 477 - DIAGNOSTIC STUART'S UREA BROTH BASE (UREA BROTH BASE) / TMV 477 - DIAGNOSTIC STUART'S UREA BROTH BASE (UREA BROTH BASE) (VEG.) / TM 1851 - UREA BROTH BASE (CHRISTENSEN) / TM 394 - UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE) / TMV 394 - UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE) (VEG.) / TM 1113 - UREA AGAR BASE (CHRISTENSEN) / TM 2410 - UREA AGAR BASE (CHRISTENSEN) (IS:5887 (Part-I)-1976) / TM 781 - MIU MEDIUM BASE / TMV 781 - MIU MEDIUM BASE (VEG.) / TM 1196 - HEMMES MEDIUM or 25ml in 975 ml TM 1554 - KOHN TWO TUBE MEDIUM NO. 1 / TM 925 - YERSINIA IDENTIFICATION BROTH BASE. Mix well and dispense in sterile tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>MedWet GmbH Südstraße 10 48163 Aachen, Germany</small>	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**

Revision: 2 March., 2022