

TS 131 – TICARCILLIN SUPPLEMENT

INTENDED USE

For selective enrichment of *Yersinia enterocolitica*.

COMPOSITION

Ingredients	Concentration
Ticarcillin	1.00mg

(per vial sufficient for 1000 ml medium)

Applied in Medium













TM 1357 - TTC BROTH BASE (ITC BROTH BASE) (IRGASAN TICARCILLIN AND POTASSIUM CHLORATE BROTH BASE)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 988 ml of sterile, molten cooled (45-50°C) TM 1357 - TTC BROTH BASE (ITC BROTH BASE) (IRGASAN TICARCILLIN AND POTASSIUM CHLORATE BROTH BASE) containing 10 ml Potassium Chlorate Supplement. Mix well and dispense into sterile tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only
Revision: 05 December., 2023