

## TS 170 – PYR REAGENT (10 ml/vl)

### INTENDED USE

For detection of Pyrrolidonyl arylamidase enzyme activity of  $\beta$ -hemolytic Streptococci.

### COMPOSITION













Ingredients	Concentration
N,N -Dimethylaminocinnamaldehyde	1.0gm
Hydrochloric acid (concentrated)	1.0ml
Distilled water	99.0ml

### INSTRUCTION FOR USE

- Inoculate test culture on PYR Agar and incubate at 35-37°C for 18-24 hours.
- Add 1 drop of PYR Reagent directly to suspected colony.
- Observe for colour change after 2 minutes.
- Formation of red colour indicates PYR positive organism.

### STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>MaxMer GmbH Max-Mer-Str. 10 48153 Ahaus, Germany</small>	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

\*For Lab Use Only  
Revision: 11 March., 2022