

## PP 005 - DISPOSABLE PETRI PLATES

### INTENDED USE

Disposable Petri plate recommended for laboratory use. For in vitro use only.

### PRODUCT FEATURES

1. Optically Clear (Non-Gamma Irradiated) for Easy Reading.
2. Packaged sterile.
3. These Disposable Polystyrene Petri Dishes are supplied in 90 mm Diameter x 15 mm height

### APPLICATION

Bioassays, Antimicrobial susceptibility tests, General lab purposes and Observation of plant germination.

### PACKAGING:

20X30 pack size.

### DISPOSAL

User must ensure safe disposal by autoclaving at a temperature of 121°C (15 psi) for 15 minutes of used medium (bacteriological materials e.g. culture media with pathogens or spores) in petri plates / unused preparations of this product on completion of work to avoid any contagion)



Quantity



Lot / Batch Number



Temperature Unit



Manufacturer



Best Before



Certification of  
Good Manufacturing Practices

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

**\*For Lab Use Only**

