

TS 160 – LANDER AF SUPPLEMENT B

INTENDED USE

For selective enrichment and transport of *Campylobacter* spp.

COMPOSITION

Ingredients	Concentration
Cyclohexamide	100mg/L

(Per vial sufficient for 1000 ml medium)

Applied in Medium


TM 1430 - LANDER TRANSPORT & STORAGE MEDIUM

INSTRUCTION FOR USE

Rehydrate the contents aseptically of 1 vial into to 1000 ml of sterile, molten TM 1430 - LANDER TRANSPORT & STORAGE MEDIUM along with Lander VPT supplement A. Mix carefully and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 11 March., 2022