

TS 308 – LACTOBACILLI SUPPLEMENT

INTENDED USE

Recommended for the selective isolation of Lactobacilli from wine.

COMPOSITION

Ingredients	Concentration
Cycloheximide	50mg

(Per vial sufficient for 500 ml medium)

Applied in Medium










TM 2379 - TOMATO JUICE MEDIUM BASE

INSTRUCTION FOR USE

Rehydrate the content of 1 vial aseptically with 2-5 ml sterile boiling distilled water. Mix well and aseptically add to 500 ml of sterile, molten cooled (45-50°C) TM 2379 - TOMATO JUICE MEDIUM BASE. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 7 March., 2022