

## TS 314 – KANAMYCIN AND VANCOMYCIN SUPPLEMENT

### INTENDED USE

An Antibiotic Supplement, recommended for the selective isolation of gram negative anaerobes.

### COMPOSITION

Ingredients	Concentration
Kanamycin	10.00 mg
Vancomycin	7.50 mg

### Applied in Medium

TM 2445 - LAKED BLOOD WITH KANAMYCIN AND VANCOMYCIN AGAR (LKV)













TM 291- SCHAEGLER AGAR

### INSTRUCTION FOR USE

One vial sufficient for 950ml medium. Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 950 ml of sterile, molten, cooled (45-50°C) Schaedler Agar along with 50 ml sterile defibrinated sheep blood and rehydrated contents of 1 vial of Vitamin K1 Supplement. Mix and pour into sterile petri plates.

### STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>MedNet GmbH Borkstrasse 10, 49163 Moenster, Germany</small>	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

**\*For Lab Use Only**  
Revision: 23 Nov., 2022