

TS 267 – ACRIFLAVIN-CEFSULODIN-VANCOMYCIN SUPPLEMENT (ACV SUPPLEMENT)

INTENDED USE

For the isolation of Enterohemorrhagic *E. coli* (EHEC).

COMPOSITION

Ingredients	Concentration
Acriflavin	50mg
Cefsulodin	40mg
Vancomycin	50mg

(Per vial sufficient for 500 ml medium)

APPLIED IN MEDIUM


TM 2020 - BUFFERED PEPTONE WATER W/ PYRUVATE

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, cooled (45-50°C) TM 2020 - Buffered peptone water w/ pyruvate. Mix well and dispense as desired.

STORAGE:

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only
Revision: 02 Feb., 2022