

TS 280 – VANCOMYCIN SUPPLEMENT

INTENDED USE

For isolation of *Cronobacter sakazakii*.

COMPOSITION

Ingredients	Concentration
Vancomycin	10mg

(per vial sufficient for 1000 ml medium)

Applied in Medium













TM 1829 - LAURYL SULPHATE TRYPTOSE BROTH, MODIFIED (ISO 22964:2006)

INSTRUCTION FOR USE

Rehydrate the content of 1 vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add it to 1000 ml of sterile, cooled (45-50°C) TM 1829 - LAURYL SULPHATE TRYPTOSE BROTH, MODIFIED (ISO 22964:2006). Mix well and pour into sterile tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only
Revision: 2 March., 2022