

TS 299 – GROWTH SUPPLEMENT I FOR MSM

INTENDED USE

A growth supplement for the enrichment of *Salmonella* species.

COMPOSITION

Ingredients	Concentration
Vitamin B12	0.100mg
Distilled Water	1ml

(per vial sufficient for 500 ml medium)

Applied in Medium


TM 2218 - MSM BROTH BASE

INSTRUCTION FOR USE

Warm up the refrigerated contents of one vial to room temperature. Aseptically add the contents of one vial of Growth Supplement I for MSM (TS 299) to 1000 ml sterile, cooled (45-50°C) TM 2218 - MSM BROTH BASE along with Growth Supplement II for MSM (TS 300). Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only
Revision: 28 Feb., 2022